

Spyder P326A Rehab

PATHMASTER

AADL #W551

OCTOBER 1, 2024

Powerchair

Freedom to move where you want when you want. The P326A Rehab offers maneuverability and affordability, meeting the needs of the individual. Perfect for both indoor and outdoor use.

- > Pathmaster Spyder 10 inch mid-wheel powerchair
- > Full suspension
- > Off-board 4A charger with 36AH batteries
- > 16" – 17" – 18" width and depth adjustable
- > Tension adjustable back with rigidizer bar
- > Right hand swing-away joystick (standard)
- > **Unoccupied** transport tie down kit (incl.)
- > Colours: red or blue



\$5,990.00
MSRP



P326 Rehab Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb°
136KG/300LBS	34"x 24"x 40"	88KG/194LBS	2 x .7HP	8KPH/5MPH	2 x36AH	F6" M10" R6"	28KM/18M	50CM/20"	6°/10%

All specifications subject to change without notice.

Order Options

- P326A
- RS326 (seat)

Colour choice

- Red
- Blue

Left hand swing-away joystick

- \$150.00 MSRP

Swing-away leg rests

- \$240.00 MSRP



www.eclipsemedical.com

#092024EN

Email: Sales@Eclipsemedical.com

Telephone: 604.540.0114

Fax: 604.540.0172