

Spyder XL P327 Rehab

PATHMASTER

AADL #W552

OCTOBER 1, 2024

Powerchair

Freedom to move where you want when you want. The P327 Rehab offers all the same benefits as its cousin the 326A, but with increased power and stability. Perfect for indoor and outdoor use.

- > Pathmaster Spyder XL 12 inch mid-wheel powerchair
- > Full suspension
- > Off-board 5A charger with 55AH batteries
- > 18" – 19" – 20" width and depth adjustable
- > Tension adjustable back with rigidizer bar
- > Right hand swing-away joystick (standard)
- > **Unoccupied** transport tie down kit (incl.)
- > Colours: red or blue



\$6,690⁰⁰
MSRP



P327 Rehab Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
205KG/450LBS	35"x 27"x 39"	110KG/242LBS	2 x .8HP	8KPH/5MPH	2 x 55AH	F6" M12" R6"	32KM/20M	53CM/21"	8°/14%

All specifications subject to change without notice.

Order Options

- P327
 RS327 (seat)

Colour choice

- Red
 Blue

Left hand swing-away joystick

- \$150.00 MSRP

Swing-away leg rests

- \$240.00 MSRP



www.eclipsemedical.com

#092024EN

Email: Sales@Eclipsemedical.com

Telephone: 604.540.0114

Fax: 604.540.0172