

VGP VOYAGER Wheelchair (Manual Wheelchair)



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____

FRAME COLOUR *Select one below* **1**

- VOYAGER PLUS FRAME - Black 1450

SEAT WIDTH *Select one below* **2**

- SW14"-20" Specify _____ N/C

SW14"-20" Light Weight Aluminum Frame

- SW21"-28" Specify _____ N/C

SW21"-28" Must Include Plus Package 1, 2 Or 3

SEAT DEPTH *Select one below* **3**

- SDF18" (By Frame) N/C
- SD16"-18" Specify _____ (By Upholstery) N/C
- SDF20" (By Frame) 350
- SD18"-20" Specify _____ (By Upholstery) N/C
- SDF22" (By Frame) 350
- SD20"-22" Specify _____ (By Upholstery) N/C

BACK *Select one below* **4**

- BH15"-20" Specify _____ N/C
- BH21"+ Specify _____ 120
- AB Angle Adj. Back N/C
- ENDO* Endomorph Canes **W ()** N/C

Specify Width

Endomorph Canes require TA or TAS Arm Rest

STROLLER HANDLE *Select one below* **5**

- SHS Removable Stroller Handle N/C
- SHR Removable Raised Stroller Handle N/C

FOOT REST *Select one below* **6**

- F6 60 Degree Hanger N/C
- F7 70 Degree Hanger N/C
- F8 80 Degree Hanger N/C
- EA Elevating Legrest **(A2)** 404
- AL Articulating Legrest **(A2)** 450

FOOT REST AND FOOT PLATES *Select one below* **7**

- PPL Plastic N/C
- PAF Angle Adjustable - Extra Large 150
- PXL Extra Large - Alum (F6A,F7A,F8A) 150
- 1FB One Piece Foot Board 80

ARM REST *Select one below* **8**

- FA Flip Back Ht. Adj./ Removable 240
- TA* T-Arm Ht. Adj./ Removable 240
- FAS Flip Back (Short) Ht. Adj./ Remove 240
- TAS* T-Arm (Short) Ht. Adj./ Removable 240
- GWAPS Gel Waterfall Arm Pads **(A1)** 150

14" - 20" Wide Light Weight Aluminum Frame Only

Please Note : A1 = Clients Responsibility

Please Note : A2 = Clinical Justification

SPECIFICATIONS

- LIGHT WEIGHT ALUMINUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- BLACK PADDED ARM PADS
- DOUBLE CROSS BRACE

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	SL 12.25	20"	5"	XL SHORT
<input type="checkbox"/>	SL 13.00"	20"	5"	XL SHORT
<input type="checkbox"/>	LL 13.75"	20"	5"	XL SHORT
<input type="checkbox"/>	L 14.50"	20"	6"	SHORT
<input type="checkbox"/>	L 14.50"	22"	5"	SHORT
<input type="checkbox"/>	L14.50"	22"	6"	SHORT
<input type="checkbox"/>	LH 15.25"	20"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	20"	6"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	6"	SHT/LONG
<input type="checkbox"/>	SH 16.00"	22"	6"	LONG

(STANDARD SEAT SETTING)

<input type="checkbox"/>	MH 16.75	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	8"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	24"	8"	LONG
<input type="checkbox"/>	TH18.25"	24"	8"	LONG
<input type="checkbox"/>	TT 19.00"*	24"	8"	LONG

***TT 19.00' requires taller caster journal**

REAR WHEELS		Select one below	9
<input type="checkbox"/>	RWT0	20" MAG Solid Treaded Tire	N/C
<input type="checkbox"/>	RWT2	22" MAG Solid Treaded Tire	N/C
<input type="checkbox"/>	RWT4	24" MAG Solid Treaded Tire	N/C
<input type="checkbox"/>	RWP2	22" Pyramid Tire	80
<input type="checkbox"/>	RWP4	24" Pyramid Tire	80
<input type="checkbox"/>	RWA2	22" Air Tire	120
<input type="checkbox"/>	RWA4	24" Air Tire	120
<input type="checkbox"/>	AQR	Axles – Quick Release (349 Lbs)	N/C
<input type="checkbox"/>	PCHR	Plastic Coated Hand Rims	N/C

14" - 20" Wide Light Weight Aluminium Frame Only

CASTER		Select one below	10
<input type="checkbox"/>	FC1A	5" X 1.25" Caster	N/C
<input type="checkbox"/>	FC1B	6" X 1.25" Caster	N/C
<input type="checkbox"/>	FC1C	8" X 1.25" Caster	N/C
<input type="checkbox"/>	FC2A	5" X 1.75" Caster	80
<input type="checkbox"/>	FC2B	6" X 1.75" Caster	80
<input type="checkbox"/>	FC2C	8" X 1.75" Caster	80
<input type="checkbox"/>	FC3A	5" X 1.25" Air Caster	80
<input type="checkbox"/>	FC3B	6" X 1.25" Air Caster	80
<input type="checkbox"/>	FC3C	8" X 1.25" Air Caster	80

WHEEL LOCK		Select one below	11
<input type="checkbox"/>	ALPL	Aluminium Pull To Lock	N/C
<input type="checkbox"/>	ALPH	Aluminium Push To Lock	N/C
<input type="checkbox"/>	STPL	Steel Pull To Lock	N/C
<input type="checkbox"/>	STPH	Steel Push To Lock	N/C
<input type="checkbox"/>	BRXA	6" Brake Extension	40
<input type="checkbox"/>	BRXB	9" Brake Extension	60

SEATING - BACKS		Select one below	12
		W () x H ()	
<input type="checkbox"/>	TAB	Tension Adjustable Back	180

ACCESSORIES			14
<input type="checkbox"/>	SBA	Seat Belt (Auto Style)	22
<input type="checkbox"/>	SBB	Seat Belt (Velcro Style)	15.75
<input type="checkbox"/>	CNC	Crutch And Cane Holder (A1)	60
<input type="checkbox"/>	IVA	Adjustable I V Pole (A1)	130
<input type="checkbox"/>	NLSP	Nylon Leg Strap Padded (A1)	50
<input type="checkbox"/>	VG PAMPL	Amputee Support Left	45
<input type="checkbox"/>	VG PAMPR	Amputee Support Right	45
<input type="checkbox"/>	VG POXY	Oxygen Tank Holder	130
<input type="checkbox"/>	VG PATR	Anti Tippers	88

Please Note : A1 = Clients Responsibility

PLUS PACKAGES		Select one below	15
---------------	--	------------------	----

- PLUS PACKAGE 1: AADL CODE: W101**
CLIENT 250 TO 349 LBS - \$880.00
INCLUDES THE FOLLOWING:
STEEL FRAME, DOUBLE CROSS BRACES WITH HEAVY DUTY BACK AND SEAT UPHOLSTERY
SEAT WIDTH INCREASE FROM 16" TO 24"
SEAT DEPTH INCREASE FROM 16" TO 22"
5,6 OR 8 X 1 3/4" HD CASTERS
*** WELDED CASTER JOURNAL ***
- PLUS PACKAGE 2: AADL CODE: W102**
CLIENT 350 TO 449 LBS - \$1140.00
INCLUDES THE FOLLOWING:
STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY
SEAT WIDTH INCREASE FROM 18" TO 26"
SEAT DEPTH INCREASE FROM 16" TO 22"
5,6 OR 8 X 1 3/4" HD CASTERS
HD STEEL FORKS
*** WELDED CASTER JOURNAL ***
- PLUS PACKAGE 3: AADL CODE: W103**
CLIENT 450 TO 600 LBS - \$1300.00
INCLUDES THE FOLLOWING:
STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY
SEAT WIDTH INCREASE FROM 20" TO 28"
SEAT DEPTH INCREASE FROM 16" TO 22"
5,6 OR 8 X 1 3/4" HD CASTERS
HD STEEL FORKS
*** WELDED CASTER JOURNAL ***
- FIXED REINFORCED ONE PIECE BACK POST **\$260.00**
(SPECIFY ANGLE _____)
- HEAVY DUTY TENSION ADJ. BACK UPH. **\$220.00**

Please Note : Reference Page 12 of the AADL Wheelchair APL for options and catalogue numbers for A wheelchairs.

Notes: