CeclipseTM Scooter & Powerchair Catalogue

SHOPRIDER

MALLMASTER

PATHMASTER

TRAILMASTER

Discovering your world again—in style

NTERCEPTOR



Powerchairs

Freedom to move where you want when you want. Manoeuvrability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.



Axis UL-8WSLA

- > Shoprider ultra-light powerchair
- > Easy as 1-2-3 to disassemble
- > 12_{AH} batteries
- > Off-board 2A charger
- > Colour: red







- > Shoprider 8 inch mid-wheel powerchair
- > Swing away joystick
- > Off-board 5A charger with 21AH batteries
- > Colours: red or blue





Navigator^M P424^M



- > Full suspension
- > Off-board 5A charger with 50AH batteries
- > Swing-away joystick
- > Colours: red, blue, black or champagne with interchangeable panels included





●●●● Navigator^L P424^L

- * Power-tilt and rehab seat available
- > Shoprider 14 inch mid-wheel powerchair
- > Full suspension
- > Off-board 5A charger with 50AH batteries
- Swing-away joystick
- Colours: red, blue, black or champagne with interchangeable panels included



		Lx Wx H							Turning	
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	Safe climb ^o
UL8wsla	113kg/250lbs	30"x 21.5"x 32"	47кg/104LBS	2 х .5нр	6 крн/ 3.9 мрн	2 х 12ан	ғ6" к9"	14км/9м	40см/15.5"	8°/14%
888wnls	113кg/250ьвѕ	37"x 22.5"x 38"	64kg/141LBS	2 х .5нр	7 крн/ 4.4 мрн	2 х 21ан	ғ6" м8" r5"	24км/15м	43см/17"	9°/15%
Р424м	136кg/300ьв	40"x 22.5"x 44"	100kg/220LBS	2 х .8нр	8крн/5мрн	2 х 50ан	ғ6" м10" R5"	34км/21м	49см/18"	9°/15%
P424L	158kg/350lbs	41"x 25"x 46"	115kg/253lbs	2х1.34нр	9.6крн/6мрн	2 х 50ан	ғ6" м14" R5"	34км/21м	57см/22.5"	9°/15%

PATHMASTER





Spyder P326A ••



- Pathmaster Spyder 10 inch mid-wheel powerchair
- Full suspension
- Off-board 4a charger with 36ah batteries
- Swing-away joystick
- > Colours: red or blue







- > Pathmaster Spyder-R 10 inch mid-wheel powerchair
- > Full suspension
- > Off-board 4A charger with 36AH batteries
- > 16" 17" 18" width and depth adjustable
- Tension adjustable back
- Swing-away joystick
- > Colours: red or blue






Spyder XL P327 ••



- Pathmaster Spyder XL 12 inch mid-wheel powerchair
- Full suspension
- Off-board 5A charger with 55AH batteries
- Swing-away joystick
- > Colours: red or blue







- Full suspension
- Off-board 5A charger with 55AH batteries
- 18" 19" 20" width and depth adjustable
- Tension adjustable back
- > Swing-away joystick
- > Colours: red or blue





		Lx Wx H							Turning	
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	Safe climb ^o
P326a	136kg/300LBS	34" x 24" x 48"	84kg/185LBS	2 х .7нр	8крн/5мрн	2 х 36ан	F6"M10"R6"	28км/18м	50см/20"	6°/10%
P326ar	136kg/300LBS	34" x 24" x 40"	88kg/194LBS	2 х .7нр	8крн/5мрн	2 х 36ан	ғ6"м10"R6"	28км/18м	50см/20"	6°/10%
P327	205kg/450LBS	35" x 27" x 50"	106кg/233LBS	2 х .8нр	8крн/5мрн	2 х 55ан	F6" M12" R6"	32км/20м	53см/21"	8°/14%
P327 _R	205kg/450lbs	35" x 27" x 39"	110kg/242lbs	2 х .8нр	8крн/5мрн	2 х 55ан	ғ6" м12" r6"	32км/20м	53см/21"	8°/14%



Compact/portable

Setting the standard in portable mobility with the world's first complete connector-less chassis design for ease of assembly & disassembly. It's as easy as 1-2-3.



Chameleon GK-83



- > Shoprider 3/4-wheel scooter
- > Lightweight and portable
- > Easy as 1-2-3 to disassemble
- > 12AH batteries with 2A charger
- > 8 inch flat free tires
- Colours: red, blue or champagne with interchangeable panels included

\$1,995.00 MSRP



Chameleon GK-84



- > Shoprider 4-wheel scooter
- > Lightweight and portable
- > Easy as 1-2-3 to disassemble
- > 12AH batteries with 2A charger
- > 8 inch flat free tires
- > Colours: red, blue or champagne with interchangeable panels included





Spirit GK10



- > Shoprider 4-wheel scooter
- > Lightweight and portable
- > Easy as 1-2-3 to disassemble
- > 35AH batteries (can be charged out of scooter) with 5A charger
- > 10 inch flat free tires
- > Nickel plated bumpers
- > Front and rear carry handles
- > Deluxe contour seat and light package
- > Colours: red or monarch blue

\$3,295.00 MSRP

		Lx Wx H							Turning	
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	Safe climbo
GK8 ³	114kg/250LBS	38"x 21.5"x 35"	41kg/89LBS	0.5нр	6 крн/ 3.8 мрн	2 х 12ан	₽7" R8"	14км/9м	84см/33"	6°/10%
GK8 ⁴	114kg/250LBS	39"x 21.5"x 35"	43kg/94LBS	0.5нр	6 крн/ 3. 8мрн	2 х 12ан	₽7" R8"	14км/9м	94см/37"	6°/10%
GK10	136kg/300lbs	48"x 22"x 38"	80kg/176LBS	0.8нр	8крн/5мрн	2 х 35ан	F10" R10"	33км/21м	133см/52"	9°/15%

MALLMASTER



Sprint S237

- > Mallmaster Sprint 3-wheel scooter
- > Lightweight and portable
- > 22AH detachable battery pack with 1.8A charger
- > Flat free tires
- > Deluxe swivel seat with sliders
- > Colours: red, blue or white







Sprint S247 •••



- > Mallmaster Sprint 4-wheel scooter
- > Lightweight and portable
- > 22AH detachable battery pack with 1.8A charger
- > Flat free tires
- > Deluxe swivel seat with sliders
- > Colours: red, blue or white

		Lx Wx H							Turning	Safe
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	climbo
S237	136кg/300ьвѕ	40"x 21"x 32"	47kg/104LBS	.6нр	7.2крн/4.5мрн	2 х 22ан	F8" R9"	17км/10.6м	86см/34"	6°/10%
S247	136кg/300ьвѕ	40"x 21"x 32"	49кg/108LBS	.6нр	7.2крн/4.5мрн	2 х 22ан	f8" r9"	17км/10.6м	120см/47"	6°/10%

Mid-size



Our classic line of scooters with the perfect blend of style and performance. All of these mid-size scooters provide the same great features with the option of 3-wheel manoeuvrability, 4-wheel stability or a combination of both with the 778EL model.

Voyager 778s



- > Shoprider 3-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > Full light package
- > Front and rear suspension
- > 36AH batteries with 5A off-board charger
- > Colours: red, silver or monarch blue





Cobra 778EL



- Shoprider 3/4-wheel scooter
- C-type delta tiller
- > Electric seat lift and full light package
- > 36AH batteries with 5A charger
- > 4-wheel stability with 3-wheel manoeuvrability
- > Colours: red, silver or monarch blue

Explorer 888SLN

- > Shoprider 4-wheel scooter
- > Delta tiller
- > Full light package
- > Full suspension, anchor hooks
- > 50AH batteries with 5A charger
- > Convenient battery access panel
- > Colours: red, black or monarch blue



\$ 3,89	95.	OO MSRP
/ -		

		Lx Wx H							Turning	
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	Safe climb ^o
778s	136kg/300LBS	45"x 22.5"x 42"	88kg/194LBS	0.87нр	8крн/5мрн	2 х 36ан	F10" R11.5"	32км/20м	125см/49"	8°/14%
778EL	136kg/300LBS	43.5"x 22"x 41"	73kg/160LBS	0.84нр	8крн/5мрн	2 х 36ан	F8" R10"	32км/20м	110см/43"	8°/14%
888sln	136kg/300lbs	47.5"x 22.5"x 42"	99kg/218LBS	0.87нр	8крн/5мрн	2 х 50ан	F10" R10"	40км/25м	146см/58"	8°/14%

TRAILMASTER



• Yeti S746

- > Trailmaster Yeti 4-wheel scooter
- > Digital LCD dashboard
- > Full light package
- > Rear mono shock suspension
- > 50ah batteries with 5a charger
- > Colours: red, blue or white







- > Trailmaster Pioneer 4-wheel scooter
- > Digital LCD dashboard
- > Dual headlights with active brake lights
- > Rear suspension
- > 50ah batteries with 5a charger
- > Colours: red, blue or silver



		Lx Wx H							Turning	
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	Safe climb ^o
S746	136kg/300LBS	48"x 24"x 48"	94kg/207LBS	1.6нр	10крн/6.25мрн	2 х 50ан	10"	28км/18м	150см/59"	10°/16%
S148	204kg/450lbs	51"x 25"x 45"	110kg/242LBS	2.3нр	12.8крн/8мрн	2 х 50ан	12"	40км/25м	148см/57"	10°/16%

Full-size



Shoprider's full-size series provides feature rich luxury in personal mobility. Experience the outdoors in unprecedented style, comfort and performance. The full-size series are the most popular scooters in the

Shoprider Canada portfolio.



TrailBlazer 889^{SL}



- > Shoprider 4-wheel scooter
- > Delta tiller with full light package
- > 1.2HP motor with 12.5 inch tires
- > Full suspension, anchor hooks
- > 50an batteries with 5a charger
- > Colours: red, blue or silver



••• TrailBlazer SE 889SL/SE



- > Special Edition Shoprider 4-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > 2HP motor with 13.5 inch tires
- > Full suspension, full light package
- > 50ah batteries with 5a charger
- > Colours: red, black or silver





LandCruiser 889XLSBN



- > Shoprider 4-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > 2HP motor with 12.5 inch tires
- > Full suspension, full light package
- > 75AH batteries with 8A charger
- > Colours: red or silver

¢_	495	00
しゅつ	495	MSRP
	, , , , , ,	MSRP

		Lx Wx H							Turning	
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	Safe climb ^o
889sl	136кg/300ьвѕ	49"x 24.5"x 51"	108kg/238lbs	1.2нр	12крн/7.45мрн	2 х 50ан	4 x 12.5"	40км/25м	150см/59"	9°/15%
889sl/se	136kg/300LBS	49"x 24.5"x 51"	112kg/246LBS	2.0нр	13.5крн/8.4мрн	2 х 50ан	4 x 13.5"	40км/25м	150см/59"	9°/15%
889xlsbn	227kg/500lbs	56"x 28"x 50"	147kg/323lbs	2.0нр	13.5крн/8.4мрн	2 х 75ан	4 x 12.5"	48км/30м	178см/70"	12°/20%

TRAILMASTER



Bigfoot 8836

- > Trailmaster Bigfoot 3-wheel scooter
- > Digital LCD dashboard
- > Full light package
- > Front & rear suspension with 400 lbs. weight capacity
- > 75AH batteries with 6A charger
- > Colours: red, blue or orange







Bigfoot S846



- > Trailmaster Bigfoot 4-wheel scooter
- > Digital LCD dashboard
- > Full light package
- > Front & rear suspension w/ 400 lbs. weight capacity
- > 75AH batteries with 6A charger
- > Colours: red, blue or orange

\$5,495.00

		Lx Wx H							Turning	Safe
Specs	Weight capacity	Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	radius	climbo
S836	182kg/400LBS	56"x 26"x 48"	127kg/280LBS	2.3нр	15крн/9.5мрн	2 х 75ан	F12" R12"	50км/32м	130см/51"	10°/16%
S846	182kg/400LBS	56"x 26"x 49"	130kg/286LBS	2.3нр	15крн/9.5мрн	2 х 75ан	F12" R12"	50км/32м	161см/63"	10°/16%

Interceptor

Premium

Interceptor S840GT •••




- > New Digital LCD Dashboard Display featuring:
 - Speedometer, Odometer, Temperature, Time and **Battery Gauges**
- > Modern L.E.D. Lighting Package includes:
 - Dual Head, Tail, Signal, Hazard and Brake lighting system
- > Robust drivetrain powered by a P&G 120 amp S - Drive controller with a high RPM four-pole motor.
- Delivers exceptional performance.
- > Infinite position Delta tiller. Fits two rear view mirrors and comes equipped with a glove box and front basket for extra storage.
- > Colours: red, blue or silver

* Batteries not included.



Interceptor S940GTX





- Speedometer, Odometer, Temperature, Time and Battery Gauges
- > Modern L.E.D. Lighting Package includes:
 - Dual Head, Tail, Signal, Hazard and Brake lighting system
- > Robust drivetrain powered by a P&G 140 amp S Drive controller with a high RPM four-pole motor. Delivers exceptional performance.
- > Delta tiller (wrap around handle bar). Fits two rear view mirrors. Comes equipped with an extra large glove box
- > Colours: red, blue or silver




* Batteries not included.



Specs	Weight capacity	Lx Wx H Dimensions	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climbo
S840 _{GT}	160kg/350lbs	51" x 27" x 49"	118kg/260LBS	2.3нр	15крн/9.5мрн	FITS 50-75	12"	30км/20м	167см/66"	10°/16%
S940 _{GTX}	204kg/450LBS	61" x 29" x 47"	165kg/363LBS	3нр	15крн/9.5мрн	FITS 75-100	F14" R16"	40км/25м	180см/71"	10°/16%


Accessories






Designs by NYACC Designs and Eclipse Medical


Eclipse Medical is the exclusive distributor of genuine NYACC accessories.



Canopy



Wscreen



Shopbasket

Rain and sun shade NYACC DESIGNS LTD.

Windscreen for delta tiller NYACC DESIGNS LTD.

Front mounted HD metal shopping basket.

Fits Explorer 888SLN only.

NYACC DESIGNS LTD.



Cane 1



Cane 2



Oxyholder


Single cane holder NYACC DESIGNS LTD.

Dual cane holder NYACC DESIGNS LTD.


Oxygen tank holder NYACC DESIGNS LTD.



Walkholder



Ramps



Enclosure

Walker holder NYACC DESIGNS LTD. Scooter and Powerchair aluminum ramps
5 ft & 6 ft single fold
8 ft multi fold
ECLIPSE MEDICAL

Full weather enclosure RRS840 (Interceptor S840GT only)



Head office:
Unit 110 - 725 Eaton Way, Delta
British Columbia, Canada V3M 6S5
Tel 604-540-0114 Fax 604-540-0172
www.eclipsemedical.com

Montreal office: 106 Nordic Avenue, Point Claire Quebec, Canada H9R 3Y2 Tel 514-428-5996 Fax 514-694-1299

