





P326A
SPYDER



P326A *SPYDER*

Features

- > Eclipse 10-inch mid-wheel powerchair delivers outstanding maneuverability in tight, compact spaces.
- > Full front and rear suspension for a smooth and comfortable ride both indoors and outdoors.
- > Dual in-line motors for enhanced efficiency, torque, range and performance.
- > Off-board 4a charger with 36AH batteries.
- > Swing-away joystick with height and width adjustable padded armrests.
- > Deluxe contoured 18" reclining high back seat with headrest.
- > Available in red or blue.  



P326A Specs

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motors	*High Speed	Battery	Tire	Charger	*Driving dist.	Turning radius	Safe climb ^o
136KG/300LBS	34"x 24"x 48"	88KG/194LBS	2 x .8HP	8KPH/5MPH	2 x 36AH	F6" M10" R6"	4AMP/HR	28KM/18M	508CM/20"	6°/10.5%

All specifications subject to change without notice.

*Speed and/or distance is variable with user weight, payload, terrain type, battery size (AH), battery charge, battery condition, ambient temperature and tire air pressure.

- Due to manufacturing tolerances and product advancements, specifications can be subject to a variance.

Dealer stamp



Discovering your world again—in style

www.eclipsemmedical.com