

PATIENT HANDLING SLINGS

made local, made right,
made to measure



Shipped in 48hrs*



The best quality, readily available, at the lowest possible cost

We are committed to being the best at what we do – supplying world-class durable medical products and superior customer service to Canada’s Health Care industry. With over eighty years cumulative experience in the moving and handling industry we understand what works and what does not.

Our designs improve on the best, most innovative features currently offered. HighStar slings fit better, last longer and are easier to use.



DELUXE
HAMMOCK




DELUXE
RECLINE



HAMMOCK



HYGIENE



POSITIONING



UNIVERSAL

 **Designed and Made in Canada by Canadians**

Manufacturing in Mississauga, Ontario allows us to fulfill orders quickly. With HighStar you are dealing with a reliable manufacturer who offers superior customer support while meeting your specific business or client needs.



HighStar Healthcare Inc.

1060 Stacey Court, Unit 1, Mississauga ON L4W 2X8

highstarhealthcare.com

HighStar's family of slings



Deluxe Hammock
Designed to provide maximum support to the upper thigh and lower back area. It can achieve neutral position, minimizing hip rotation. Includes best features from universal and hammock slings.



Deluxe Recline
A full body sling designed to eliminate pressure points, preventing folds and creases. Excellent for clients with amputations and poor skin integrity.




Hammock
Designed to provide maximum support to the upper thigh and lower back area. It can achieve neutral position, minimizing hip rotation.



Hygiene
Designed for toileting with maximum support. Self-adjusting waist strap allows precise and comfortable positioning during a transfer.




Positioning Designed for supine transfers and repositioning of a client in bed.



Universal
General transfer sling that provides excellent trunk and thigh support. Great for clients with limited upper body tone.

EasyMove Straps

Accessory to any loop style sling with head support. They can be used as extensions to existing straps or as a counter balance against an actively moving client



Adjustable lanyards

All our head support slings can be ordered with the regular loop style straps or these adjustable lanyards.



Sling Size Guide

Size	Height	Weight
Child	< 4'	< 55 lbs
Junior	< 4'	< 110 lbs
Small*	4' – 5'6"	95-150 lbs
Medium*	5' – 6'	150-240 lbs
Large*	5'6" – 7'	240-350 lbs
X-Large	5'6" – 7'	350-450 lbs
XX-Large	> 6.5'	> 400 lbs

* Shipped in 48 hours for orders under 5 units of the same SKU

Safe Working Load Chart

Model	Material	Size	SWL
Hammock & Universal	Woven, Quilted & Mesh	Ch/Jr S/M/L XL/XXL	450 lbs. 650 lbs. 1000 lbs.
	Spacer	Ch/Jr S/M/L XL/XXL	450 lbs. 650 lbs. 800 lbs.
Deluxe Hammock & Deluxe Recline	Spacer	Ch/Jr S/M/L XL/XXL	450 lbs. 650 lbs. 800 lbs.
Hygiene	Woven	Ch/Jr S/M/L XL/XXL	450 lbs. 650 lbs. 1000 lbs.
Positioning	Spacer & Woven	One Size	1000 lbs.



HighStar
HEALTHCARE INC.

HighStar Healthcare Inc.

1060 Stacey Court, Unit 1, Mississauga ON L4W 2X8

highstarhealthcare.com