

eclipse™ medical inc.

PATHMASTER ***SPYDER XL P327***



- **12" Mid-Wheel Powerchair**
- **Adjustable Suspension**
- **Swing Away Joystick**
- **20 x 18 Deluxe High Back Seat**
- **55 ah Batteries**
- **Colours: Red & Blue**

		P327
Dimensions L x W x H	cm (in)	90 (35.4) x 69 (27.1) x 127 (50)
Total Weight w/Batteries	kg (lb.)	106 (233)
Battery	_V_Ah x Pcs	12V / 55 AH x 2
Charger	Off Board	5 amp
Tires	Front	2 x 6"
	Middle	2 x 12"
	Rear	2 x 6"
Motor	DC24V	575W (Max) x 2
Max. Speed	km/hr. (mph)	8 (5)
Driving Distance*	km (mile)	32 (20)
Weight Capacity	kg (lb.)	204 (450)
Suspension		Yes (Adjustable)
Controller		Dynamic Shark 50A

Specifications subject to change without notice

*Maximum driving distance is based on a new fully charged battery depending on weight, terrain and temperature.

Discovering your world again—in style