

# **eclipse™** *PATHMASTER* **SPYDER XL P327R**



- 12" Mid-Wheel Powerchair
- Adjustable Suspension
- Swing Away Joystick
- 18"-19"-20" Width & Depth Adj
- Tension Adjustable Back
- 55 ah Batteries
- Colours: Red & Blue

		<b>P327R</b>
<b>Dimensions L x W x H</b>	<b>cm (in)</b>	<b>90 (35) x 69 (27) x 99 (39)</b>
<b>Total Weight w/Batteries</b>	<b>kg (lb.)</b>	<b>110 (242)</b>
<b>Battery</b>	<b>_V_Ah x Pcs</b>	<b>12V / 55 AH x 2</b>
<b>Charger</b>	<b>Off Board</b>	<b>5 amp</b>
<b>Tires</b>	<b>Front</b>	<b>2 x 6"</b>
	<b>Middle</b>	<b>2 x 12"</b>
	<b>Rear</b>	<b>2 x 6"</b>
<b>Motor</b>	<b>HP / DC24V</b>	<b>.8 / 575W (Max) x 2</b>
<b>Max. Speed</b>	<b>km/hr. (mph)</b>	<b>8 (5)</b>
<b>Driving Distance*</b>	<b>km (mile)</b>	<b>32 (20)</b>
<b>Weight Capacity</b>	<b>kg (lb.)</b>	<b>205 (450)</b>
<b>Suspension</b>		<b>Yes (Adjustable)</b>
<b>Controller</b>		<b>Dynamic Shark 50A</b>

Specifications subject to change without notice

Ver 1.2

\*Maximum driving distance is based on a new fully charged battery depending on weight, terrain and temperature.

Discovering your world again—in style